

# 3M Food Safety Environmental Listeria Monitoring Solution



## Zero In On *Listeria*



3M™ Petrifilm™  
Environmental Listeria Plate

3M™ Quick Swab




Monitor environmental *Listeria* levels with consistent, frequent testing to identify hot spots in your plant and track your progress. Cost effectively test more often in more locations. Zero in on the source of *Listeria* with a complete solution from 3M Food Safety.


- Quantitative results allow you to identify hot spots in the plant
- Rapid results confirmed in 27–31 hours from sample collection
- Safer method—no enrichment, no culture transfers—less risk of cross-contamination
- No specialised equipment needed
- Easy-to-use solution to make faster decisions, track progress against *Listeria* over time, and react to problem areas




## Fast, Easy and More Informative




1. Sample using swab or sponge.



2. Add buffered peptone water. Wait one hour.



3. Inoculate and incubate for 26–30 hours.



4. Read and interpret.

## End-to-end solution that works together for consistent, reliable results

### Product Solutions

#### Microbial Testing—For detection of *Listeria* in environmental samples

Product	Catalog Number	Quantity
3M™ Petrifilm™ Environmental Listeria Plate	6447	50 plates
	6448	200 plates

#### Sample Collection

Product	Catalog Number	Quantity
3M™ Quick Swab	6432	50 swabs
	6433	250 swabs
3M™ Swab-Sampler w/4mL Neutralising Buffer	RS9604NB	100 swabs
3M™ Swab-Sampler w/10mL Neutralising Buffer	RS96010NB	100 swabs
3M™ Enviro Swab	ENVS WB25	25 swabs
3M™ Sponge-Stick, 1.5" x 3" Sponge with 10mL Neutralising Buffer	SSL10NB	100 sponges (5 bags of 20)
3M™ Hydrated-Sponge, 1.5" x 3" Sponge with 10mL Neutralising Buffer	HS10NB	100 sponges (5 bags of 20)

#### Plating

Product	Catalog Number	Quantity
3M™ Petrifilm™ Spreader for Environmental Listeria Plate	6498	2 each

For more information visit  
[www.3M.com/FoodSafety](http://www.3M.com/FoodSafety) or give us a call

**3M**

**ThermoFisher**  
 SCIENTIFIC

For customer service, call 1300-735-292  
 To fax an order, use 1800-067-639  
 Visit us online: [www.thermofisher.com.au](http://www.thermofisher.com.au)

©2011 Thermo Fisher Scientific Inc. All rights reserved. A.B.N. 52 058 390 917



3M and Petrifilm are trademarks of 3M.